

# Amigo Arc

Keeps your child in the loop



child  
friendly  
hearing  
care




## Let children engage

At school, obstacles such as noise, distance and reverberation can make concentrating very difficult. By sending the teacher's voice directly to the child's ears, FM receivers help them use their energy not just on keeping up, but on staying one step ahead of the game.

### What is Amigo Arc?

Amigo Arc is a universal neckloop receiver that acts as a bridge between an FM transmitter and the child's hearing aid. Amigo Arc works with any instrument with a telecoil - including cochlear implants and bone anchored devices. It can also be used with headphones by children with listening or attention difficulties, to help them hear better and concentrate more easily.



FM > Arc > telecoil

FM > Arc > headphones



## Get more out of music

Amigo Arc offers easy connection to mp3 players - with a bass boost that makes music sound richer than ever. No other FM neckloop can offer this!

When it's time to watch TV, children will get far more out of their favourite programmes, as the extra FM boost makes the sound quality much clearer.



## Small, light - and safe

To the untrained eye, Amigo Arc looks much like an mp3 player. Children can show it off in one of five trendy colours - from cool white to chilli red - or conceal it completely under a shirt or sweater. Either way, it's extremely safe to wear, thanks to a safety release that automatically opens if the neckloop should get caught on something.



Personalise Amigo Arc with colourful shells that are easy to click on and off.



**child  
friendly  
hearing  
care**

It takes a truly dedicated approach to help children with hearing problems achieve their full potential. That's why we deliver the solutions and services that professionals and caregivers need to provide children the opportunities they deserve. This is what child-friendly hearing care is all about.



100557/03.10

[www.amigofm.com](http://www.amigofm.com)  
[www.oticonchildren.com](http://www.oticonchildren.com)

**oticon**  
PAEDIATRICS